

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 505

Burettes, 50 ml, lat. valve., Schellb., AS, conf. cert.

Batch: 18-26




Serial:

Measuring Nr. 79513

Batch size:

31

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	$\pm 0,0500$ ml	$\pm 0,0500$ ml	$\pm 0,0500$ ml
Average value:	49,9840 ml	29,9858 ml	10,0104 ml
Deviation:	0,0115 ml	0,0127 ml	0,0147 ml
Precision (CV):	0,0230 %	0,0424 %	0,1471 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

06.08.2026

Balance:

Sartorius -
LE244S-OCE

DE-M