

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 014 253

Burettes, 25 ml, str. stopc., amber glass, AS, conf. cert.

Batch: 43-25




Serial:

Measuring Nr. 79033

Batch size:

60

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	25,0030 ml	15,0031 ml	4,9969 ml
Deviation:	0,0014 ml	0,0025 ml	0,0022 ml
Precision (CV):	0,0055 %	0,0168 %	0,0440 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

15.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M