

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 122 505

Pellet-burettes, 50 ml, lat. needle valve, AS, conf. cert.

Batch: 18-26




Serial:

Measuring Nr. 79465

Batch size:

25

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0173 ml	30,0368 ml	10,0268 ml
Deviation:	0,0126 ml	0,0059 ml	0,0024 ml
Precision (CV):	0,0252 %	0,0198 %	0,0239 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

06.08.2026

Balance:

Sartorius -
LE244S-OCE

DE-M