

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 124 108

Pellet-burettes, 10 ml, amber glass, AS, conf. cert.




Batch: 23-21

Serial:

Measuring Nr. 72081

Batch size: 10

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	$\pm 0,0200$ ml	$\pm 0,0200$ ml	$\pm 0,0200$ ml
Average value:	9,9893 ml	5,9902 ml	1,9994 ml
Deviation:	0,0027 ml	0,0033 ml	0,0043 ml
Precision (CV):	0,0272 %	0,0551 %	0,2156 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.03.2025

Balance:

Sartorius -
LE244S-OCE

DE-M