

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 383 025

Micro burettes, 2 ml, lat. needle valve, AS, conf. cert.

Batch: 40-25




Serial:

Measuring Nr. 78940

Batch size:

10

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	$\pm 0,0100$ ml	$\pm 0,0100$ ml	$\pm 0,0100$ ml
Average value:	1,9966 ml	1,1971 ml	0,4007 ml
Deviation:	0,0031 ml	0,0018 ml	0,0019 ml
Precision (CV):	0,1568 %	0,1534 %	0,4728 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

08.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M