

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 001 BL

Grad. cylinders, 10:0,2 ml, blue grad, A, conf. cert.




Batch: 08-25

Serial:

Measuring Nr. 72179

Batch size: 500

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
Average value:	2,0083 ml	6,0028 ml	9,9972 ml
Deviation:	0,0123 ml	0,0123 ml	0,0111 ml
Precision (CV):	0,6160 %	0,2052 %	0,1114 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

05.03.2025

Balance:

Sartorius -
LE244S-OCE

DE-M