

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL

Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.




Batch: 26-26

Serial:

Measuring Nr. 79079

Batch size: 500

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	10,0060 ml	30,0418 ml	50,0730 ml
Deviation:	0,0537 ml	0,0604 ml	0,0752 ml
Precision (CV):	0,5366 %	0,2013 %	0,1505 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

17.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M