

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf.cert.




Batch: 19-25

Serial:

Measuring Nr. 72897

Batch size: 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	$\pm 0,2500$ ml	$\pm 0,2500$ ml	$\pm 0,2500$ ml
Average value:	9,9838 ml	30,0205 ml	50,0594 ml
Deviation:	0,0465 ml	0,0506 ml	0,0476 ml
Precision (CV):	0,4651 %	0,1688 %	0,0953 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

07.05.2025

Balance:

Sartorius -
LE244S-OCE

DE-M