

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf. cert.




Batch: 27-26

Serial:

Measuring Nr. 79280

Batch size: 280

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	$\pm 0,2500$ ml	$\pm 0,2500$ ml	$\pm 0,2500$ ml
Average value:	9,9287 ml	30,0298 ml	50,1138 ml
Deviation:	0,0548 ml	0,0407 ml	0,0439 ml
Precision (CV):	0,5485 %	0,1357 %	0,0878 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

29.07.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M