

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005

Graduated cylinders, 50:1 ml, A, conf. cert.




Batch: 27-26

Serial:

Measuring Nr. 79314

Batch size: 500

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	10,0275 ml	30,1163 ml	50,1791 ml
Deviation:	0,0561 ml	0,0461 ml	0,0397 ml
Precision (CV):	0,5612 %	0,1537 %	0,0793 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

31.07.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M