

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 13-25

Serial:

Measuring Nr. 72460

Batch size: 500

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	19,9533 ml	60,0396 ml	100,1467 ml
Deviation:	0,1211 ml	0,1063 ml	0,1128 ml
Precision (CV):	0,6055 %	0,1772 %	0,1128 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

01.04.2025

Balance:

Sartorius -
LE244S-OCE

DE-M