

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 32-26

Serial:

Measuring Nr. 79539

Batch size: 500

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	20,0579 ml	60,0787 ml	100,1333 ml
Deviation:	0,0826 ml	0,0739 ml	0,0830 ml
Precision (CV):	0,4128 %	0,1232 %	0,0830 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

11.08.2026

Balance:

Sartorius - LA
1200 S

DE-M