

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010

Graduated cylinder, 100:1,0 ml, A, conf. cert.




Batch: 13-25

Serial:

Measuring Nr. 72530

Batch size: 1000

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	$\pm 0,5000$ ml	$\pm 0,5000$ ml	$\pm 0,5000$ ml
Average value:	20,0210 ml	60,0399 ml	100,0750 ml
Deviation:	0,1300 ml	0,1015 ml	0,0881 ml
Precision (CV):	0,6502 %	0,1691 %	0,0881 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

07.04.2025

Balance:

Sartorius -
LE244S-OCE

DE-M