

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 025 USP

**Graduated cylinders, 250:2 ml, A, USP, conf. cert.**




**Batch:** 17-25

**Serial:**

**Measuring Nr.** 72798

**Batch size:** 247

|                        |                   |                    |                    |
|------------------------|-------------------|--------------------|--------------------|
| <b>Volume:</b>         | <b>50,0000 ml</b> | <b>150,0000 ml</b> | <b>250,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 1,0000$ ml   | $\pm 1,0000$ ml    | $\pm 1,0000$ ml    |
| <b>Average value:</b>  | 50,0061 ml        | 150,0068 ml        | 250,0048 ml        |
| <b>Deviation:</b>      | 0,1846 ml         | 0,2117 ml          | 0,2251 ml          |
| <b>Precision (CV):</b> | 0,3692 %          | 0,1412 %           | 0,0901 %           |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

24.04.2025

**Balance:**

Sartorius - LA  
1200 S

DE-M