

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 100 BL

Grad. cylinder, 1000 ml, NS 45, blue grad., A, conf. cert.




Batch: 28-25

Serial:

Measuring Nr. 73640

Batch size: 49

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	$\pm 5,0000$ ml	$\pm 5,0000$ ml	$\pm 5,0000$ ml
Average value:	198,5667 ml	599,8420 ml	999,9493 ml
Deviation:	0,6525 ml	0,4512 ml	0,6846 ml
Precision (CV):	0,3262 %	0,0752 %	0,0685 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

14.07.2025

Balance:

Sartorius -
MSE14202S-
OCE-D0

DE-M