

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 012 W

**Vol. flasks, 10 ml, NS 10, white grad., A, conf. cert.**

**Batch:** 18-25

**Serial:**

**Measuring Nr.** 72967

**Batch size:**

149


**Volume:** 10,0000 ml

**Tolerance:**  $\pm 0,0400$  ml

**Average value:** 10,0037 ml

**Deviation:** 0,0033 ml

**Precision (CV):** 0,0329 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

13.05.2025

**Balance:**

Sartorius -  
LE244S-OCE

DE-M