

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 020 BL

**Vol. flask, 200 ml, NS 14, blue grad., A, conf. cert.**

**Batch:** 27-26

**Serial:**

**Measuring Nr.** 79204

**Batch size:**

1191


**Volume:** 200,0000 ml

**Tolerance:**  $\pm 0,1500$  ml

**Average value:** 199,9911 ml

**Deviation:** 0,0205 ml

**Precision (CV):** 0,0102 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

23.07.2026

**Balance:**

Sartorius - LA  
1200 S

DE-M