

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 100

Volumetric flasks, 1000 ml, NS 24, A, conf. cert.

Batch: 27-26

Serial:

Measuring Nr. 78756

Batch size:

806


Volume: 1000,0000 ml

Tolerance: $\pm 0,4000$ ml

Average value: 1000,0816 ml

Deviation: 0,0538 ml

Precision (CV): 0,0054 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

01.07.2026

Balance:

Sartorius -
MSE14202S-
OCE-D0

DE-M