

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 051 BL

Grad. pipettes, 0,5:0,01 ml, blue grad., AS, conf. cert.




Batch: 23-25

Serial:

Measuring Nr. 73643

Batch size: 280

Volume:	0,5000 ml	0,3000 ml	0,1000 ml
Tolerance:	$\pm 0,0060$ ml	$\pm 0,0060$ ml	$\pm 0,0060$ ml
Average value:	0,4989 ml	0,2981 ml	0,0996 ml
Deviation:	0,0011 ml	0,0012 ml	0,0006 ml
Precision (CV):	0,2159 %	0,3847 %	0,6275 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

14.07.2025

Balance:

Sartorius -
LE244S-OCE

DE-M