

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 510 051

**Graduated pipettes, 0,5:0,01 ml, AS, conf. cert.**




**Batch:** 23-25

**Serial:**

**Measuring Nr.** 73645

**Batch size:** 280

|                        |                  |                  |                  |
|------------------------|------------------|------------------|------------------|
| <b>Volume:</b>         | <b>0,5000 ml</b> | <b>0,3000 ml</b> | <b>0,1000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0060$ ml  | $\pm 0,0060$ ml  | $\pm 0,0060$ ml  |
| <b>Average value:</b>  | 0,4995 ml        | 0,2990 ml        | 0,0997 ml        |
| <b>Deviation:</b>      | 0,0010 ml        | 0,0011 ml        | 0,0009 ml        |
| <b>Precision (CV):</b> | 0,2067 %         | 0,3748 %         | 0,8916 %         |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

15.07.2025

**Balance:**

Sartorius -  
LE244S-OCE

DE-M