

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 101

Grad. pipettes, 1:0,01 ml, AS, conf. cert.




Batch: 20-26

Serial:

Measuring Nr. 79201

Batch size: 2484

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
Tolerance:	$\pm 0,0070$ ml	$\pm 0,0070$ ml	$\pm 0,0070$ ml
Average value:	0,9994 ml	0,6020 ml	0,2010 ml
Deviation:	0,0027 ml	0,0009 ml	0,0006 ml
Precision (CV):	0,2719 %	0,1434 %	0,2904 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

22.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M