

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 202

Grad. pipettes, 2:0,02 ml, AS, conf. cert.

Batch: 23-25




Serial:

Measuring Nr. 73646

Batch size:

500

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	$\pm 0,0100$ ml	$\pm 0,0100$ ml	$\pm 0,0100$ ml
Average value:	1,9968 ml	1,1993 ml	0,4021 ml
Deviation:	0,0028 ml	0,0035 ml	0,0011 ml
Precision (CV):	0,1416 %	0,2901 %	0,2848 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

15.07.2025

Balance:

Sartorius -
LE244S-OCE

DE-M