

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 510 505 BL

Grad. pipettes, 5:0,05 ml, blue grad., AS, conf. cert.

Batch: 23-25

Serial:




-

Measuring Nr. 73658

Batch size:

280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	5,0037 ml	2,9997 ml	0,9983 ml
Deviation:	0,0027 ml	0,0021 ml	0,0031 ml
Precision (CV):	0,0539 %	0,0708 %	0,3108 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

16.07.2025

Balance:

Sartorius -
LE244S-OCE

DE-M