

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 515 101

**Grad. pipette, 1 : 0,01 ml, AS, conf. cert.**




**Batch:** 48-25

**Serial:**

**Measuring Nr.** 76029

**Batch size:** 497

Volume:	1,0000 ml	0,6000 ml	0,2000 ml
<b>Tolerance:</b>	± 0,0070 ml	± 0,0070 ml	± 0,0070 ml
<b>Average value:</b>	0,9969 ml	0,5962 ml	0,1996 ml
<b>Deviation:</b>	0,0014 ml	0,0012 ml	0,0007 ml
<b>Precision (CV):</b>	0,1410 %	0,1932 %	0,3395 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	12.12.2025
<b>Balance:</b>	Sartorius - MSE224S-0CE -DU

DE-M