

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 510

Grad. pipettes, 5:0,1 ml, AS, conf. cert.




Batch: 20-25

Serial:

Measuring Nr. 72969

Batch size: 280

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	$\pm 0,0300$ ml	$\pm 0,0300$ ml	$\pm 0,0300$ ml
Average value:	4,9856 ml	2,9901 ml	1,0007 ml
Deviation:	0,0057 ml	0,0037 ml	0,0023 ml
Precision (CV):	0,1135 %	0,1220 %	0,2302 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

14.05.2025

Balance:

Sartorius -
LE244S-OCE

DE-M