

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 515 810

Grad. pipettes, 20:0,1 ml, AS, conf. cert.




Batch: 25-26

Serial:

Measuring Nr. 78724

Batch size: 150

Volume:	20,0000 ml	12,0000 ml	4,0000 ml
Tolerance:	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
Average value:	19,9775 ml	11,9880 ml	4,0014 ml
Deviation:	0,0185 ml	0,0112 ml	0,0116 ml
Precision (CV):	0,0925 %	0,0931 %	0,2903 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 835

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

30.06.2026

Balance:

Sartorius -
LE244S-OCE

DE-M