

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 011

Pipette f.encym.analysis., 0,1:0,001 ml, AS, conf. cert.

Batch: 25-26




Serial:

Measuring Nr. 78874

Batch size:

217

Volume:	0,1000 ml	0,0600 ml	0,0200 ml
Tolerance:	$\pm 0,0030$ ml	$\pm 0,0030$ ml	$\pm 0,0030$ ml
Average value:	0,0997 ml	0,0593 ml	0,0190 ml
Deviation:	0,0006 ml	0,0009 ml	0,0007 ml
Precision (CV):	0,5783 %	1,5471 %	3,6478 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

03.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M