

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 517 202

Pipettes f. encym. analysis, 2:0,02 ml, AS, conf. cert.




Batch: 26-26

Serial:

Measuring Nr. 79034

Batch size: 143

Volume:	2,0000 ml	1,2000 ml	0,4000 ml
Tolerance:	$\pm 0,0100$ ml	$\pm 0,0100$ ml	$\pm 0,0100$ ml
Average value:	2,0012 ml	1,2012 ml	0,3990 ml
Deviation:	0,0023 ml	0,0026 ml	0,0023 ml
Precision (CV):	0,1153 %	0,2171 %	0,5664 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to DIN 12699

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

14.07.2026

Balance:

Sartorius -
LE244S-OCE

DE-M