

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 535 025

Volumetric pipettes, 2,5 ml, AS, conf. cert.

Batch: 26-26

Serial:

Measuring Nr. 79636

Batch size:

499


Volume: 2,5000 ml

Tolerance: $\pm 0,0100$ ml

Average value: 2,5016 ml

Deviation: 0,0030 ml

Precision (CV): 0,1200 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

19.08.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M