

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 535 150 USP

**Volumetric pipette, 15 ml, USP, AS, conf. cert**

**Batch:** 19-25


**Serial:**

**Measuring Nr.** 73365

**Batch size:**

144

|                        |                   |
|------------------------|-------------------|
| <b>Volume:</b>         | <b>15,0000 ml</b> |
| <b>Tolerance:</b>      | $\pm 0,0300$ ml   |
| <b>Average value:</b>  | 15,0028 ml        |
| <b>Deviation:</b>      | 0,0050 ml         |
| <b>Precision (CV):</b> | 0,0331 %          |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 648

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

23.06.2025

**Balance:**

Sartorius -  
LE244S-OCE

DE-M