

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 902

**Burettes, 100 ml, str. stopc., solid plug, AS, conf. cert.**

**Batch:** 38-16




**Serial:**

**Measuring Nr.:** 29177

**Batch size:**

1

<b>Volume:</b>	<b>100,000 ml</b>	<b>60,000 ml</b>	<b>20,000 ml</b>
<b>Tolerance:</b>	$\pm 0,1000$ ml	$\pm 0,1000$ ml	$\pm 0,1000$ ml
<b>Measured value:</b>	99,9192 ml	59,9399 ml	19,9884 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

0

**Date:**

09.04.2019

**Balance:**

Sartorius -  
LE244S-OCE

DE-M