

Certificate of performance



Article: 3 010 053

Burettes, 5 ml, straight stopcock, AS, conf. cert.




Batch: 10-20

Serial:

Measuring Nr. 76833

Batch size: 8

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	5,0027 ml	2,9954 ml	0,9995 ml
Deviation:	0,0033 ml	0,0011 ml	0,0004 ml
Precision (CV):	0,0665 %	0,0373 %	0,0424 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	02.03.2026
Balance:	Sartorius - MSE224S-0CE -DU

