

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 103 B

Burettes, 10 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 26-24

Serial:

Measuring Nr. 74981

Batch size: 60

| Volume: | 10,000 ml | 6,000 ml | 2,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0200 ml | ± 0,0200 ml | ± 0,0200 ml |
| Average value: | 10,0004 ml | 5,9996 ml | 1,9981 ml |
| Deviation: | 0,0073 ml | 0,0027 ml | 0,0033 ml |
| Precision (CV): | 0,0725 % | 0,0453 % | 0,1628 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 16.10.2025

Balance: Sartorius -
LE244S-OCE

DE-M