

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 103

Burettes, 10 ml, straight stopcock, AS, conf. cert.




Batch: 30-17

Serial:

Measuring Nr. 13409

Batch size: 50

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0099 ml	6,0016 ml	2,0016 ml
Deviation:	0,0034 ml	0,0022 ml	0,0012 ml
Precision (CV):	0,0342 %	0,0375 %	0,0581 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	23.08.2017
Balance:	Sartorius - LE244S-OCE

