

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 253 B

Burettes, 25 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 13-18

Serial:

Measuring Nr.: 20292

Batch size: 400

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	25,0009 ml	14,9798 ml	4,9815 ml
Deviation:	0,0049 ml	0,0025 ml	0,0019 ml
Precision (CV):	0,0196 %	0,0168 %	0,0384 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	03.05.2018
Balance:	Sartorius - LE244S-OCE

