

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 010 253

Burettes, 25 ml, str. stopcock, PTFE, AS, conf. cert.




Batch: 18-21

Serial:

Measuring Nr. 47063

Batch size: 60

Volume:	25,0000 ml	15,0000 ml	5,0000 ml
Tolerance:	± 0,0300 ml	± 0,0300 ml	± 0,0300 ml
Average value:	24,9998 ml	14,9948 ml	5,0012 ml
Deviation:	0,0068 ml	0,0036 ml	0,0046 ml
Precision (CV):	0,0273 %	0,0240 %	0,0919 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	01.06.2021
Balance:	Sartorius - LE244S-OCE

