

Certificate of performance



Article: 3 010 503 B

Burettes, 50 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 17-25

Serial:

Measuring Nr. 75252

Batch size: 298

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0169 ml	29,9902 ml	10,0078 ml
Deviation:	0,0134 ml	0,0135 ml	0,0137 ml
Precision (CV):	0,0268 %	0,0451 %	0,1368 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	11.11.2025
Balance:	Sartorius - LE244S-OCE

