

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 503

Burettes, 50 ml, str. stopcock, PTFE, AS, conf. cert.




Batch: 20-23

Serial:

Measuring Nr. 63562

Batch size: 97

Volume:	50,0000 ml	30,0000 ml	10,0000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0016 ml	29,9950 ml	10,0058 ml
Deviation:	0,0124 ml	0,0112 ml	0,0136 ml
Precision (CV):	0,0248 %	0,0372 %	0,1365 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	29.06.2023
Balance:	Sartorius - LE244S-OCE

DE-M