

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 010 503

**Burettes, 50 ml, str. stopcock, PTFE, AS, conf. cert.**




**Batch:** 02-24

**Serial:**

**Measuring Nr.** 70193

**Batch size:** 99

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Average value:	50,0121 ml	30,0068 ml	10,0100 ml
Deviation:	0,0158 ml	0,0086 ml	0,0098 ml
Precision (CV):	0,0315 %	0,0286 %	0,0980 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	77
Date:	13.11.2024
Balance:	Sartorius - LE244S-0CE

