

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 010 903 B

Burette, 100 ml, str. stopc., amber grad., AS, conf. cert.




Batch: 50-21

Serial:

Measuring Nr. 52099

Batch size: 60

| Volume: | 100,000 ml | 60,000 ml | 20,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,1000 ml | ± 0,1000 ml | ± 0,1000 ml |
| Average value: | 99,9943 ml | 59,9759 ml | 20,0072 ml |
| Deviation: | 0,0328 ml | 0,0214 ml | 0,0238 ml |
| Precision (CV): | 0,0328 % | 0,0356 % | 0,1188 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 10.02.2022 |
| Balance: | Sartorius - LE244S-OCE |

DE-M