

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 052

**Burettes, 5 ml, straight stopc., glass, AS, conf. cert.**




**Batch:** 07-26

**Serial:**

**Measuring Nr.** 77624

**Batch size:** 10

| Volume:         | 5,0000 ml   | 3,0000 ml   | 1,0000 ml   |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value:  | 5,0009 ml   | 3,0000 ml   | 1,0023 ml   |
| Deviation:      | 0,0035 ml   | 0,0026 ml   | 0,0026 ml   |
| Precision (CV): | 0,0704 %    | 0,0853 %    | 0,2557 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                           |
|-----------|---------------------------|
| Operator: | 7                         |
| Date:     | 23.04.2026                |
| Balance:  | Sartorius -<br>LE244S-OCE |

DE-M