

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 053

Burettes, 5 ml, str. stopc., PTFE plug, AS, conf. cert.




Batch: 24-16

Serial:

Measuring Nr. 12578

Batch size: 2

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
Tolerance:	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
Average value:	4,9995 ml	2,9964 ml	0,9955 ml
Deviation:	0,0010 ml	0,0001 ml	0,0021 ml
Precision (CV):	0,0199 %	0,0047 %	0,2058 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	02.08.2017
Balance:	Sartorius - LE244S-OCE

