

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 053

**Burettes, 5 ml, str. stopc., PTFE- plug, AS, conf. cert.**




**Batch:** 09-21

**Serial:**

**Measuring Nr.** 46192

**Batch size:** 29

| Volume:         | 5,0000 ml   | 3,0000 ml   | 1,0000 ml   |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,0100 ml | ± 0,0100 ml | ± 0,0100 ml |
| Average value:  | 5,0026 ml   | 2,9984 ml   | 1,0001 ml   |
| Deviation:      | 0,0031 ml   | 0,0021 ml   | 0,0009 ml   |
| Precision (CV): | 0,0614 %    | 0,0704 %    | 0,0917 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                                   |
|-----------|-----------------------------------|
| Operator: | 4                                 |
| Date:     | 26.04.2021                        |
| Balance:  | Sartorius -<br>MSE224S-0CE<br>-DU |

