

# Certificate of performance



**Article:** 3 012 053

**Burettes, 5 ml, str. stopc., PTFE- plug, AS, conf. cert.**




**Batch:** 40-22

**Serial:**

**Measuring Nr.** 60411

**Batch size:** 38

Volume:	5,0000 ml	3,0000 ml	1,0000 ml
<b>Tolerance:</b>	± 0,0100 ml	± 0,0100 ml	± 0,0100 ml
<b>Average value:</b>	4,9952 ml	2,9946 ml	0,9980 ml
<b>Deviation:</b>	0,0024 ml	0,0018 ml	0,0010 ml
<b>Precision (CV):</b>	0,0483 %	0,0615 %	0,0997 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	77
<b>Date:</b>	22.11.2022
<b>Balance:</b>	Sartorius - LE244S-OCE

