

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 103

**Burettes, 10 ml, str. stopc., PTFE-plug, AS, conf. cert.**




**Batch:** 36-16

**Serial:**

**Measuring Nr.** 15518

**Batch size:** 30

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
<b>Tolerance:</b>	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
<b>Average value:</b>	9,9999 ml	5,9977 ml	1,9991 ml
<b>Deviation:</b>	0,0053 ml	0,0029 ml	0,0026 ml
<b>Precision (CV):</b>	0,0529 %	0,0486 %	0,1275 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

9

**Date:**

27.10.2017

**Balance:**

Sartorius -  
LE244S-OCE

DE-M