

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 103

Burettes, 10 ml, str. stopc., PTFE-plug, AS, conf. cert.




Batch: 32-18

Serial:

Measuring Nr. 23810

Batch size: 150

Volume:	10,000 ml	6,000 ml	2,000 ml
Tolerance:	± 0,0200 ml	± 0,0200 ml	± 0,0200 ml
Average value:	10,0044 ml	5,9986 ml	1,9974 ml
Deviation:	0,0048 ml	0,0049 ml	0,0048 ml
Precision (CV):	0,0483 %	0,0818 %	0,2412 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	07.09.2018
Balance:	Sartorius - LE244S-OCE

DE-M