

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 103

Burettes, 10 ml, str. stopc., PTFE-plug, AS, conf. cert.

Batch: 45-25




Serial:

Measuring Nr. 78049

Batch size:

97

Volume:	10,0000 ml	6,0000 ml	2,0000 ml
Tolerance:	$\pm 0,0200$ ml	$\pm 0,0200$ ml	$\pm 0,0200$ ml
Average value:	9,9916 ml	5,9866 ml	2,0017 ml
Deviation:	0,0034 ml	0,0028 ml	0,0026 ml
Precision (CV):	0,0336 %	0,0466 %	0,1314 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

11.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M