

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 012 253

**Burettes, 25 ml, str. stopc., PTFE plug, AS, conf. cert.**




**Batch:** 10-26

**Serial:**

**Measuring Nr.** 77953

**Batch size:** 97

|                        |                         |                         |                         |
|------------------------|-------------------------|-------------------------|-------------------------|
| <b>Volume:</b>         | <b>25,0000 ml</b>       | <b>15,0000 ml</b>       | <b>5,0000 ml</b>        |
| <b>Tolerance:</b>      | $\pm 0,0300 \text{ ml}$ | $\pm 0,0300 \text{ ml}$ | $\pm 0,0300 \text{ ml}$ |
| <b>Average value:</b>  | 25,0067 ml              | 15,0057 ml              | 5,0065 ml               |
| <b>Deviation:</b>      | 0,0065 ml               | 0,0088 ml               | 0,0054 ml               |
| <b>Precision (CV):</b> | 0,0260 %                | 0,0590 %                | 0,1076 %                |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

05.05.2026

**Balance:**

Sartorius -  
LE244S-OCE

DE-M