

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 502 TR

Burette, 50 ml, with funnel, AS, conf. cert.

Batch: 03-17




Serial:

Measuring Nr. 14855

Batch size:

1

Volume:	50,000 ml	30,000 ml	10,000 ml
Tolerance:	± 0,0500 ml	± 0,0500 ml	± 0,0500 ml
Measured value:	50,0045 ml	30,0263 ml	10,0283 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

11.10.2017

Balance:

Sartorius -
LE244S-OCE

DE-M