

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 012 505

Burettes, 50 ml, lat. stopc., Schellbach, AS, conf. cert.




Batch: 49-16

Serial:

Measuring Nr. 16819

Batch size: 20

| Volume: | 50,000 ml | 30,000 ml | 10,000 ml |
|-----------------|-------------|-------------|-------------|
| Tolerance: | ± 0,0500 ml | ± 0,0500 ml | ± 0,0500 ml |
| Average value: | 50,0080 ml | 29,9901 ml | 9,9883 ml |
| Deviation: | 0,0162 ml | 0,0103 ml | 0,0139 ml |
| Precision (CV): | 0,0325 % | 0,0342 % | 0,1390 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 385

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|-----------|---------------------------|
| Operator: | 9 |
| Date: | 03.01.2018 |
| Balance: | Sartorius - LE244S-OCE |

DE-M